



## SINGLE WORK TUBE FURNACES



- Our VCTF range of single work tube furnaces utilises a wire wound ceramic work tube, providing excellent temperature uniformity down the length of the work tube
- Vecstar's VCTF tube furnaces are all equipped with a mesh heat shield as standard, this keeps the operator safe while the furnace is in use. Only ultra 'high grade' insulation materials are used to ensure our furnaces are as energy efficient as possible.
- All VCTF furnaces are available in either a horizontal or vertical position.

### ABOUT VECSTAR

We build furnaces, ovens and hot plates for industry, laboratories and educational institutions. During the 30+ years Vecstar have been established, a comprehensive range of advanced electric heating equipment has been developed.

### CONTACT US

W: [www.vecstar.com](http://www.vecstar.com)

E: [info@vecstar.com](mailto:info@vecstar.com)

UK Tel:

+44 (0)1246 260094

USA Tel:

1-877-842-4652 (Toll Free)

### ADDRESS

Head Office Address:

Units 11 & 12

Foxwood Road

Dunston Trading Estate

Chesterfield

S41 9RF

Model	Max Temp	Max. continuous working temperature	Chamber Diameter (mm)	Heated Length (mm)	External Dimensions (mm)			Power (kW)	Holding Power (kW)	T/C Type
					Height	Width	Depth			
VCTF1	1200°C	1150°C	15	130	320	220	150	0.4	0.16	K
VCTF2	1200°C	1150°C	25	200	550	310	250	1.25	0.5	K
VCTF3	1200°C	1150°C	38	330	550	310	380	2.25	0.9	K
VCTF4	1200°C	1150°C	50	450	550	310	500	2.5	1	K
VCTF5	1200°C	1150°C	65	600	550	310	650	3	1.2	K
VCTF6	1200°C	1150°C	75	700	550	310	750	4	1.6	K
VCTF7	1200°C	1150°C	100	950	550	310	1000	5	2	K