



BESPOKE TUBE FURNACES

ABOUT VECSTAR

We build furnaces, ovens and hot plates for industry, laboratories and educational institutions. During the 30+ years Vecstar have been established, a comprehensive range of advanced electric heating equipment has been developed.

CONTACT US

W: www.vecstar.com

E: info@vecstar.com

UK Tel:

+44 (0)1246 260094

USA Tel:

1-877-842-4652 (Toll Free)

ADDRESS

Head Office Address:

Units 11 & 12

Foxwood Road

Dunston Trading Estate

Chesterfield

S41 9RF



- Vecstars team of designers and engineers have extensive experience designing and manufacturing bespoke tube furnaces. Over the years Vecstar have created numerous bespoke tube furnace designs and solutions such as rotating tube furnaces, stainless steel furnaces; for use in clean rooms and multi-positional tube furnaces for maximum flexibility.

- A number of examples are shown below, please contact us for more information.

- All our high temperature tube furnaces are compatible with our gas tight ends, this allows inert gasses to be introduced and contained within the chamber

