



STANDARD HOT PLATES

Our industrial hot plates are robust and hard wearing. We offer hot plates with a continuous operating temperature of up to 300°C (572°F) with the thermostat option and up to 400°C (752°F) with a digital controller fitted. Our hot plates are available with either Mild Steel, Stainless Steel or Aluminium work surfaces. The body of our hot plates are entirely built from Stainless Steel, ideal for laboratories and clean room environments.

We offer 4 standard hotplate sizes sizes:

Size 1 : 100 x 275 x 8mm (3.9" x 10.8" x 0.3")

Size 2 : 305 x 305 x 8mm (12" x 12" x 0.3")

Size 3 : 305 x 610 x 8mm (12" x 24" x 0.3")

Size 4 : 305 x 915 x 8mm (12" x 36" x 0.3")

We build large and bespoke hotplates on request.

ABOUT VECSTAR

We build furnaces, ovens and hot plates for industry, laboratories and educational institutions. During the 30+ years Vecstar have been established, a comprehensive range of advanced electric heating equipment has been developed.

CONTACT US

W: www.vecstar.com

E: info@vecstar.com

UK Tel:

+44 (0)1246 260094

USA Tel:

1-877-842-4652 (Toll Free)

ADDRESS

Head Office Address:

Units 11 & 12

Foxwood Road

Dunston Trading Estate

Chesterfield

S41 9RF

